

# Quantitative Analysis of Trace Amount of Butadiene Rubber in High Impact Polystyrene (HIPS)

High impact polystyrene (HIPS) is a polystyrene (PS) copolymerized with a few percents of butadiene rubber in order to increase impact resistance of PS. Pyrolysis GC is used as a tool for quality control of HIPS. Quantitative analysis of a trace amount of butadiene present in HIPS is described here using Double-Shot Pyrolyzer®.

Fig. 1 shows a pyrogram of HIPS obtained at 550°C. Butadiene (B) and styrene (S), monomer components of HIPS, have been detected. Table 1 shows peak area ratios of B and S (B/S) in the pyrograms obtained repeatedly. The excellent reproducibility of 0.62% was obtained. In the actual quantitative analysis, various mixing ratios of samples are used to produce a calibration curve.

Table 1. Peak Area Ratios of Butadiene (B) and Styrene (S) and Reproducibility

| n       | B/S (%) |
|---------|---------|
| 1       | 1.249   |
| 2       | 1.244   |
| 3       | 1.246   |
| 4       | 1.254   |
| 5       | 1.233   |
| Average | 1.245   |
| RSD (%) | 0.62    |




Fig. 1. Pyrogram of HIPS

Pyrolysis temp: 550°C, Split ratio: 1/50, GC oven temp: 40–320°C at 20°C/min  
 Separation column: Ultra ALLOY+5 (5% diphenylpolysiloxane), Length 30m, Id 0.25mm, Film thickness 0.25µm  
 Sample amount : 100µg, Detector : Hydrogen flame detector (FID)

**Keywords :** Impact Resistance, Polystyrene, Butadiene, Quantitative Analysis, Reproducibility

**Products used :** Multi-functional pyrolyzer, UA-5

**Applications :** Plastics and Rubber Industry, General Polymer Analysis

**Related technical notes :**

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**  
**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
<http://www.frontier-lab.com/>