

Analysis of Food Wrap Films Using Double-Shot Pyrolyzer®

Part 3: Analysis of Polypropylene (PP) + Nylon by EGA GC/MS Technique

Polypropylene (PP) + Nylon food wrap films were analyzed using EGA-GC/MS technique. Fig. 1 shows an EGA profile obtained by programmed pyrolysis from 40°C~600°C at 30°C/min. Fig. 2 shows results of GC/MS analysis of temperature zones A (100~320°C), and B (200~600°C) employing MicroJet Cryo-Trap (MJT-1030E). In Zone A, volatile acetic acid, and fatty acids and their derivatives as plasticizer were found. In Zone B, olefinic hydrocarbons of C₆, C₉, C₁₂, and C₁₅ derived from pyrolysis of polypropylene, and ε-caprolactam, monomer of nylon-6, were found.




Fig. 1 EGA Profile of Polypropylene + Nylon

Pyrolysis temp: 40~600°C (30°C/min), carrier gas: He
Interface: deactivated metal capillary column (length: 2.5m, id: 0.15mm)
Injection port pressure: 50kPa




Fig. 2 Analysis Results of Zones A, B, and C of Polypropylene + Nylon

Carrier gas: He, column flow rate: 1ml/min, total carrier gas flow rate: 40ml/min, separation column: Ultra ALLOY-5 (5% diphenyl dimethyl polysiloxane), length: 30m, id: 0.25mm, film thickness: 0.25μm, GC oven temp: 40°C (1min hold) ~ 320°C (20°C/min), injection port temp: 320°C, Cryo trap temp: -196°C, sample: 0.25cm²

Reference: Hosaka et.al., 49th Japan Analytical Society Meeting (2000)

Keyword : Food Wrap Film, Evolved Gas, Plasticizer

Applications : Film manufacturer, Food producer, General polymer analysis

Please forward your inquiries via our web page at :
<http://www.frontier-lab.com/>, or send us a fax message.

R&D and manufactured by :

Frontier Laboratories Ltd.

Phone: (81)24-935-5100 Fax: (81)24-935-5102

Double-Shot Pyrolyzer® is a registered trademark of Frontier Laboratories Ltd.

Your dealer: