

Analysis of a plastic bag containing calcium carbonate

Part 2: Pyrolysis (Py)-GC/MS and thermal desorption (TD)-GC/MS

[Background] In the previous note (PYA3-040E), evolved gas analysis (EGA)-MS of a shopping bag made of calcium carbonate and thermoplastics was reported (Fig. 1). In this note, the qualitative analysis of polymer types by pyrolysis (Py)-GC/MS and qualitative and quantitative analysis of additives by thermal desorption (TD)-GC/MS were conducted.

[Experimental] A Py-GC/MS system with a Multi-Shot Pyrolyzer (EGA/PY-3030D) directly interfaced with the GC inlet was used for the measurements. A freeze-ground shopping bag was used as a sample and placed in an Eco-Cup LF. For pyrolysis, the sample was introduced into the pyrolyzer furnace heated at 600 °C for flash pyrolysis to obtain a pyrogram. For thermal desorption, volatiles were cryo-trapped using a MicroJet Cryo-Trap, and the trapped compounds were separated and analyzed to obtain a TD chromatogram. The additive of Irgafos 168 was determined by the standard addition method.

[Results] The pyrogram obtained from the flash pyrolysis of the sample is shown in Fig. 2 (a). The profile of the pyrogram indicates that the main component is polyethylene. TD chromatogram shown in Fig. 2 (b) indicates the presence of additives such as stearate metal salt, antioxidant, and polyethylene wax. Assuming that Irgafos168 undergoes the same reaction in TD-GC/MS as its oxide does, the quantification of Irgafos168 was conducted by the standard addition method, and the quantified value was 196 ppm.




Fig. 1 EGA thermogram of sample





Fig. 2 (a) Pyrogram and (b) TD chromatogram of a sample

Furnace temp.: 600 °C, Thermal desorption temp.: 100 - 400 °C (20 °C/min), GC injector temp.: 300 °C
 GC oven temp.: 40 (2 min hold) - 320 °C (20 °C/min), Split ratio: 1/100 (pyrolysis), 1/20 (thermal desorption)
 Separation column: UA⁺-5 (5 % diphenyl 95 % dimethylpolysiloxane), L=30 m, i.d.=0.25 mm, df=0.25 μm
 Column flow rate: 1.0 mL/min (He), MS scan range: *m/z* 29 - 800, Sample amount: 0.193 mg (pyrolysis), 1.108 mg (thermal desorption)

Keywords : Py-GC/MS, TD-GC/MS, Calcium carbonate, Limestone, Shopping bag

Products used : Multi-functional pyrolyzer, Eco-Cup LF, UA⁺-5, Vent-free GC/MS adapter, F-Search, MicroJet Cryo-Trap

Applications : General polymer analysis

Related technical notes : [PYA3-040E](#), [PYA3-028E](#), [PYA3-029E](#), [PYA1-119E](#), [PYA1-120E](#), [PYA1-121E](#)

Please forward your inquiries via our web page or send us a fax message.

R&D and manufactured by :
Frontier Laboratories Ltd.

Phone: (81)24-935-5100 Fax: (81)24-935-5102
www.frontier-lab.com