

## Trace analysis of acrylic copolymers using pyrolysis (Py)-GC/MS with F-Splitless injection Part 1: Comparison of pyrograms obtained with various injection methods

**[Background]** Pyrolysis (Py)-GC/MS usually uses a split injection method with a minimum sample amount of about 1 µg. In recent years, the demand for trace analysis has led to the need for the analysis of a more smaller sample amount. In response to this demand, we have developed the F-Splitless injection method<sup>1)</sup>, which enables highly sensitive detection of pyrolyzates by splitless injection using a Multi-Functional Splitless Sampler (MFS). In this report, as a basic study on the trace analysis of acrylic copolymers by Py-GC/MS using the F-Splitless injection method, pyrograms obtained by various injection methods were compared.

**[Methods]** A specified amount of a dichloromethane solution of acrylic copolymers consisting of seven monomers (Table 1) was put in a sample cup and the solvent was allowed to evaporate. Then, measurements were done using a Py-GC/MS system shown in Fig. 1.

**[Results]** The pyrograms obtained by different injection methods are shown in Fig. 2. In the split injection method, peaks of all seven monomers were detected, while in the splitless injection method, the pyrogram profile was very different from that obtained by the split injection method. This was due to the secondary reactions occurring under a slow carrier gas flow. On the other hand, the F-Splitless injection method afforded a pyrogram that was similar to that obtained by the split injection.

This is because in the F-Splitless injection method, the total flow rate is increased to about 30 mL/min by operating the MFS during pyrolysis, thereby suppressing the secondary reactions. In the next note (PYA1-161E), the high sensitivity of the F-Splitless injection method is described.






Fig. 1 Py-GC/MS system

(A) Split injection (1/50, 50 µg; equivalent to 1 µg introduced into separation column)



(B) Splitless injection (1 µg)



(C) F-Splitless injection (1 µg)




Fig. 2 Pyrograms of acrylic copolymers obtained by different injection methods.

Furnace temp.: 600 °C, Injector press.: 150 kPa, Pre-column: UA<sup>+</sup>-5 (5 % diphenyl 95 % dimethylpolysiloxane, L=2 m, i.d.=0.25 mm, df=0.25 µm), Separation column: UA<sup>+</sup>-5 (5 % diphenyl 95 % dimethylpolysiloxane, L=30 m, i.d.=0.25 mm, df=0.25 µm), GC oven: 40 °C (3 min hold) - 320 °C (20 °C/min, 3 min hold), MS scan range: m/z 29 - 550, MS scan rate: ca. 3 scan/s, Sample amount/injection method: see main text.

1) K. Tei et al., *J. Anal. Appl. Pyrolysis*, 2022, 168, 105707.

**Keywords** : Acrylic resin, Trace analysis, F-Splitless Injection, High sensitivity analysis, Flash Pyrolysis (Py)-GC/MS

**Products used** : Multi-Shot Pyrolyzer, Multi-Functional Splitless Sampler, MicroJet Cryo-Trap, Auto-Shot Sampler, UA<sup>+</sup>-5, Eco-Cup LF, F-Search, Vent-free GC/MS adapter

**Applications** : General polymer analysis, Quality assurance, Material analysis, Foreign materials analysis

**Related technical notes** : [PYA1-161E \(Part 2\)](#), [PYA1-162E \(Part 3\)](#), [PYT-037E](#), [PYA1-154E](#)

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**  
**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
[www.frontier-lab.com](http://www.frontier-lab.com)