

Analysis of Polybutylene Terephthalate (PBT) by Reactive Pyrolysis

When analyzing condensation polymers such as polybutylene terephthalate (PBT) by Py-GC technique, reactive pyrolysis in the presence of tetramethyl ammonium hydroxide (TMAH) gives constituent monomer of methyl ester.

Fig. 1 shows a pyrogram obtained by flash pyrolysis of PBT, and Fig 2 shows a pyrogram obtained by reactive pyrolysis in the presence of TMAH. Flash pyrolysis technique gave products arising from decomposition and decarboxylation of ester group, but no monomer. On the other hand, reactive pyrolysis gave PBT constituent monomer of dimethyl derivatives of terephthalic acid and mono and dimethyl derivatives of 1,4-butanediol.




Fig. 1 Pyrogram Obtained by Flash Pyrolysis of PBT




Fig. 2 Pyrogram Obtained by Reactive Pyrolysis of PBT

Analytical conditions: carrier gas: He, Injection port pressure: 103kPa, Split ratio: 1/60, Separation column: Ultra ALLOY⁺-5 (5% diphenyldimethylpolysiloxane) Length: 30m, Id: 0.25mm, Film thickness: 0.25μm, GC oven temp: 38°C~300°C (20°C/min), GC injection port temp: 320°C
Material excerpted from "5. A few recent applications of Py-GC", Kiura, Wakabayashi (Mitsubishi Rayon), 2nd Pyrolysis Gas Chromatography Seminar (hosted by Frontier Lab Ltd.)

Keywords : Polybutylene terephthalate, PBT, Reactive Pyrolysis, Tetramethyl ammonium hydroxide, TMAH

Products used : Multi-functional pyrolyzer, UA-5

Applications : General Polymer Analysis

Related technical notes :

Please forward your inquiries via our web page or send us a fax message.

R&D and manufactured by :
Frontier Laboratories Ltd.

Phone: (81)24-935-5100 Fax: (81)24-935-5102
<http://www.frontier-lab.com/>

©: A registered trademark of Frontier Laboratories Ltd.