


# Analysis Examples Using Carrier Gas Switcher

## Part 3: Pyrolysis of Polycarbonate (PC) in Air

Using Carrier Gas Switcher (CGS-1050E) , Selective Sampler (SS-1010E), and MicroJet Cryo-Trap (MJT-1030E), flash pyrolysis of polycarbonate (PC) was performed at 550°C both in air and He atmosphere. Pyrolyzates were analyzed by GC/MS. Fig. 1 compares pyrograms obtained both in air and He atmosphere. In the pyrogram obtained in He atmosphere, a large amount of bis-phenol A, monomer of PC, was found, in addition to phenols such as phenol and p-cresol. On the other hand, in the pyrogram obtained in air, a large amount of carbon dioxide was observed with a trace amount of phenols. The results demonstrates that when PC was pyrolyzed, a major portion of it was decomposed and oxidized in air at high temperature.



**Fig. 1 Comparison of Pyrograms of PC Obtained in Air and He Atmosphere**

Pyrolysis temp: 550°C, Carrier gas: Helium, Column flow rate: 1 ml/min, Carrier gas total flow rate: 60 ml/min  
 Separation column: Ultra ALLOY-5 (5% diphenyldimethylpolysiloxane), Length: 30 m, Id: 0.25 mm, Film thickness: 0.25 µm  
 GC oven temperature: 40°C (1min hold) ~ 320°C (20 °C/min), Injection port temp: 320°C, Sample size: 3.0 µg

Reference: Hosaka et al. 5th Polymer Analysis Symposium, II-4, p43-44 (2000)

**Keywords :** Polycarbonate, Carrier Gas Switcher, Pyrolysis in air

**Products used :** Multi-functional pyrolyzer, Carrier Gas Selector, Selective Sampler, MicroJet Cryo-Trap, UA-5

**Applications :** General Polymer Analysis, Environmental analysis

**Related technical notes :**

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**  
**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
<http://www.frontier-lab.com/>

©: A registered trademark of Frontier Laboratories Ltd.