




# Conversion of ethanol to ethylene using Tandem μ-Reactor GC/MS system – Part 2

**[Background]** Part 1 of this series of Technical Notes (RXA-001E) used the online catalytic conversion of ethanol to ethylene to illustrate how the Tandem μ-Reactor GC/MS can be used to study, real time, catalytic conversion processes. This note describes the conversion of ethanol to ethylene using a series of specific temperature steps. The reaction products formed at each temperature were analyzed by GC/MS.

**[Experimental]** The system configuration is the same as that used in the real time study described in RXA-001E - see Fig. 1. Ethanol vapors in the carrier gas continuously flow through the 1<sup>st</sup> Reactor (100°C isothermal). The 2<sup>nd</sup> Reactor was heated in discrete steps: 100°C, 250°C, 300°C and 400°C (Fig. 2). The reaction products present at each temperature were directed to the separation column using a Selective Sampler (SS-1010E)\*<sup>1</sup>. The reaction products are trapped at the head of the GC column prior to starting the oven temperature program using a MicroJet Cryo-Trap (MJT-1035E)\*<sup>2</sup>. The zeolite catalyst, H-ZSN-5 (3φ x 10 mm) is packed in a quartz tube located in 2<sup>nd</sup> reactor.

**[Results]** The chromatograms of the catalytic reaction products at each temperature are shown in Fig. 3. At 100°C, ethanol the major component; however, as the temperature increases, the conversion to ethylene increases. Plots of the peak areas of products vs. temperature are shown in Fig. 4. As the temperature increases, the level of diethyl ether (DE) formed from the dehydration of ethanol increases, and reaches a maximum at 250°C. At temperatures higher than 250°C the concentrations of water and ethylene increase due to the dehydration of DE. This simple stepwise catalytic conversion illustrates one of the many operating modes of the Tandem μ-Reactor GC/MS system.



\*1), 2) Please see the catalogue of each accessory for details.

Ref.: R. Freeman et al., *J. Anal. Appl. Pyrolysis*, 111 (2015) 41-46

**Keywords :** Catalyst screening, Evaluation, Conversion of ethanol to ethylene

**Products used :** Tandem μ-Reactor, μ-Reactor, Selective Sampler, MicroJet Cryo-Trap, Vent-free GC/MS adapter, UA-1

**Applications :** Catalyst screening and evaluation

**Related technical notes :** RXT-001E, RXA-001E

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**

**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
<http://www.frontier-lab.com/>