


# Analysis of Polywax 655

## Using Ultra ALLOY®-SIM Column for Simulated Distillation

Ultra ALLOY®-SIM (stationary phase : Dimethylpolysiloxane) column has been designed for simulated distillation analysis.

This is a low bleeding nonpolar column, especially useful for simulated distillation analysis, in which the high thermal stability of the Ultra ALLOY® columns (see Ultra ALLOY® Column Tech Note UAI-003E) is best utilized. Fused silica capillary columns are always vulnerable to column breakage at high temperatures due to thermal degradation of polyimide outer wall coating in the air; however, Ultra ALLOY®-SIM metal column is free from such problems.

Specifications: Ultra ALLOY®-SIM (Dimethylpolysiloxane), 5m (0.53mmi.d.) 0.1µm, Temp.range: - 60 ~430°C



**Fig.1 Analysis of Polywax 655**

Column: Ultra ALLOY®-SIM (Dimethylpolysiloxane), 5 m (0.5 mm i.d.), 0.1 µm  
 Oven temp.: 40 ~ 20 °C/min ~ 430°C(1 min), Injector: Split 1/20 at 360°C  
 Detector: FID at 430°C, Carrier gas: He, 0.5 psi  
 Sample size: 1.0%, 1.0 µl (solvent : cyclohexane)

**Keywords :** Polywax 655, High Temperature Analysis

**Products used :** Multi-functional pyrolyzer, Ultra ALLOY®-SIM

**Applications :** Analysis of High Boiling Components, Crude oil / Lubricant Analysis

**Related technical notes :**

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**  
**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
<http://www.frontier-lab.com/>