



Micro-puncher

Convenient Tool for Sampling Film-Like Materials

Simply apply the tip of the tool against a thin or film-like material, and rotate the tool gently, then push out a sample piece.

■ Variations

id (mm) 0.5 0.75 1.25 2.0 3.0 4.0 5.0



Pushing the button of Micro-puncher, and stainless steel (or a white PTFE) nozzle pushes out a sample piece.

■ Sampling of papers



■ Sampling of leaves



■ Micro-punchers

| Product name | Part number | Specifications | | |
|----------------------|-------------|--|--------|--------|
| | | Contents | id(mm) | Color |
| Micro-puncher 050S | FMP-0.50S | puncher (3 pcs), Cutting Mat L (1 sheet) | 0.5 | White |
| Micro-puncher 075S | FMP-0.75S | puncher (3 pcs), Cutting Mat L (1 sheet) | 0.75 | White |
| Micro-puncher 125S | FMP-1.25S | puncher (3 pcs), Cutting Mat L (1 sheet) | 1.25 | White |
| Micro-puncher 200S | FMP-2.00S | puncher (3 pcs), Cutting Mat L (1 sheet) | 2.0 | Gray |
| Micro-puncher 300S | FMP-3.00S | puncher (3 pcs), Cutting Mat L (1 sheet) | 3.0 | White |
| Micro-puncher 400S | HMP-4.00S | puncher (3 pcs), Cutting Mat L (1 sheet) | 4.0 | Brown |
| Micro-puncher 500S | HMP-5.00S | puncher (3 pcs), Cutting Mat L (1 sheet) | 5.0 | Orange |
| Micro-puncher 03 Kit | FMP-03KIT | 3 pcs each of 0.75 / 1.25 / 2.0 (mm id) punchers and 2 sheets of Cutting Mat L | | |
| Micro-puncher 04 Kit | HMP-04KIT | 3 pcs each of 3.0 / 4.0 / 5.0 (mm id) punchers and 2 sheets of Cutting Mat L | | |

■ Cutting Mat

| Product name | Part number | Specifications | | |
|---------------|-------------|---------------------|------------|-------------|
| | | Quantity in package | D x W (mm) | Color |
| Cutting Mat L | FCM-6577 | 5 | 150 x 220 | Light green |

- There is no kit that contains Micro-puncher with 0.5 mm id.
- Sets of 10 Micro-punchers with a same id are available.



FRONTIER LABORATORIES LTD.

<http://www.frontier-lab.com/>

4-16-20, Saikon, Koriyama, Fukushima, 963-8862 Japan Phone: 81-24-935-5100 Fax: 81-24-935-5102